

EDGC NICE service

Non-Invasive Prenatal Testing



EONE-DIAGNOMICS
Genome Center






Objectives:

- What is NIPT or NIPS?
- History of NIPT
- Recommendations of NIPT
- EDGC methodology
- Case examples
- Argument for all chromosome testing

Maternal Age Related Risk for Chromosome Abnormalities

Maternal Age	First Trimester		Second Trimester		Live-birth	
	Down Syndrome	All	Down Syndrome	All	Down Syndrome	All
20	1 in 1152		1 in 1211		1 in 1477	
21	1 in 1125		1 in 1184			
22	1 in 1110		1 in 1168			
23	1 in 1090		1 in 1147			
24	1 in 1064		1 in 1120			
25	1 in 1032		1 in 1085			
26	1 in 978		1 in 1029			
27	1 in 928		1 in 997			
28	1 in 856		1 in 901			
29	1 in 775		1 in 827			
30	1 in 868		1 in 733			
31	1 in 591		1 in 632			
32	1 in 494		1 in 536			
33	1 in 401		1 in 435			
34	1 in 315	1 in 346				
35	1 in 240	1 in 114	1 in 265	1 in 141	1 in 353	1 in 179
36	1 in 179	1 in 87	1 in 197	1 in 111	1 in 267	1 in 149
37	1 in 131	1 in 66	1 in 147	1 in 88	1 in 199	1 in 123
38	1 in 96	1 in 51	1 in 108	1 in 70	1 in 148	1 in 105
39	1 in 71	1 in 38	1 in 80	1 in 56	1 in 111	1 in 81
40	1 in 53	1 in 28	1 in 60	1 in 44	1 in 85	1 in 63
41	1 in 41	1 in 22	1 in 47	1 in 35	1 in 67	1 in 49
42	1 in 32	1 in 17	1 in 38	1 in 28	1 in 54	1 in 39
43	1 in 27	1 in 13	1 in 31	1 in 22	1 in 45	1 in 39
44	1 in 22	1 in 10	1 in 26	1 in 18	1 in 39	1 in 31
45	1 in 19	1 in 8	1 in 23	1 in 14	1 in 35	1 in 24







ACMG recommends:





ACMG recommends:

- “Laboratories should provide readily visible and clearly stated detection rate (DR), clinical specificity (SPEC), positive predictive value (PPV), and negative predictive value (NPV) for conditions being screened, in pretest marketing materials, and when reporting laboratory results to assist patients and providers in making decisions and interpreting results.
- Laboratories should not offer to screen for Patau, Edwards, and Down syndromes if they cannot report DR, SPEC, and PPV for these conditions
- All laboratories should include a clearly visible fetal fraction on NIPS reports.
- All laboratories should establish and monitor analytical and clinical validity for the fetal fraction.
- All laboratories should specify the reason for a no-call when reporting NIPS results



	Screening	How to	Since when	How long	Detection Rate(%)
NIPT	NICE	Non-Invasive	From 10 weeks	7~10 days	>99%
Conventional Blood Test	Triple Screen	Non-Invasive	From 11-13 weeks	2 days	67~71%
	Quadruple Screen		From 11-13 weeks		79~81%
Integrated Screening Test	Integrated Screen	Non-Invasive	From 11-13 weeks From 11-13 weeks	4~5 weeks	94~96%
Cell Culture Test	Chorionic Screen	Invasive	From 11-13 weeks	1~2 weeks	>99%
	Amniocentesis		From 11-13 weeks		







Table1. Differences between NICE® and targeted sequencing methods

NICE® - Massively Parallel Sequencing	Targeted Sequencing
All chromosomal abnormalities can be detected	Only major chromosomal abnormalities can be detected
Unlike other MPS-based NIPT tests, it reports using 21 z-score thresholds	Reported as a risk score similar to serum screening
There is also no effect due to differences between ethnicities	Depending on the SNP, it may be affected by differences between ethnicities
Amplification of fetal-derived cfDNA/maternal-derived cfDNA by size selection method using paired-end sequencing	Inability to isolate fetal-derived cfDNA and maternal-derived cfDNA

Sensitivity (False-positive rates)	EDGC NICE®	Sequenom MaterniT21	Illumina Verifi®	Ariosa Harmony	Natera Panorama
Trisomy 21 Down syndrome	>99% (<0.01%)	99.1% (0.2%)	>99.9% (0.3%)	>99% (<0.1%)	99% (0%)
Trisomy 18 Edwards syndrome	96.5% (<0.01%)	96.9% (<0.01%)	97.3% (0.4%)	96.7% (<0.1%)	94.1% (<0.1%)
Trisomy 13 Patau syndrome	92.31% (<0.01%)	89.3% (0.3%)	87.5% (0.1%)	80% (0.05%)	>99% (0%)
Monosomy X Turner syndrome	>99.99% (<0.01%)	94.7% (0.5%)	95% (1.0%)	96.7% (unreported)	94.7% (<0.1%)
Sex chromosome Trisomies	>99.99% (0%)	>99.9%	67-100%	67-100%	73.1%
Female	>99% (0%)	97.9% (0.5%)	97.6% (0.8%)	>99% (unreported)	>99.9% (0%)
Male	>99% (0%)	99.4% (2.1%)	99.1% (1.1%)	>99% (<0.01%)	>99.9% (0%)

Double Z-score



Multi-Z

Is the core technology to detect fetus chromosomal abnormalities









Figure 1. Normalization between chromosomes. Normalized value between two chromosomes is calculated by dividing the value of interested chromosome by that of each chromosome.




Size Selection







XO Case : Not detected (Maternal Mosaicism)






Trisomy 18 : 2-step confirmation







Multi-Z 10mb heatmap




Stair-Matrix
(in-house algorithm)




WiseCondor



	Full	Partial	SUM
T7	47	1	48
T19	23	0	23
T16	12	1	13
T8	10	3	13
T20	9	1	10
T10	7	2	9
T15	7	1	8
T3	8	0	8
T9	8	0	8
T14	5	2	7
T22	7	0	7
T11	5	1	6
T2	3	0	3
T17	2	0	2
T6	2	0	2
T1	0	1	1
T12	1	0	1
T4	1	0	1
T5	1	0	1
SUM	158	13	171





21,644 samples





UPD	Parent	Syndrome/Disorder	Phenotype
6	Paternal	Transient neonatal diabetes mellitus	IUGR, neonatal diabetes
7	Maternal	Russell-Silver	IUGR/FTT, dysmorphic
11	Paternal	Beckwith-Wiedemann	Omphalocele, organomegaly, neonatal hypoglycemia, Wilms tumor
11	Maternal	Russell-Silver	IUGR/FTT, dysmorphic
14	Paternal	Temple syndrome	IUGR, dysmorphic
14	Maternal	Kagami-Ogata syndrome	Bell-shaped thorax, developmental retardation, dwarfisms, dysmorphic
15	Maternal	Prader-Willi	Obesity, dysmorphic, ID
15	Paternal	Angelman	ID, dysmorphic
20	Maternal	Growth failure, hyperactivity	IUGR/FTT
20	Paternal	Pseudohypoparathyroidism 1b	Pseudohypoparathyroidism



T7 partial, FF: 9.75 10w 4d





T9, FF: 5.5 17w 6d
47,X(),+9[14]/46,X()[16]




T16 FF: 14.98 17w 3d
47,X(),+16[3]/46.XO[37]




T22, FF: 10.92 17w 0d
47,X(),+22





- 31 yo female
- NICE at 11w5d - negative, FF 10.3%
- Repeated NICE at 24w1d- negative, FF 21.2% (no nasal bone)




- No amnio- baby born 47, XY+21
- No Nasal Bone ultrasound
- Repeat NICE should not have been performed
- CPM cause



Sample	Original	Fetus	Mom	Final Sex	Fetal Fraction	FFGap	Sample Bias QC	Size QC	Final result	MD results	Pregnancy type
ON221102255.LP22120005	Not Detected	XO Suspected	Not Detected	XX	8.29	NA	PASS	PASS	XO Detected	Not Detected	Singleton



- Mom and fetus Xp21.1 - 11.3



Test Option

	NICE® LITE	NICE® BASIC	NICE® PREMIUM
T13, T18, T21	✓	✓	✓
T9, T16, T22		✓	✓
All Chromosome			✓
*8 Microdeletions			
*116 Microdeletions			
*Sex Chromosome Disorder			

* Any or all can be added to LITE, BASIC, or PREMIUM service

NICE® Test Report



Sample Information

Sample Type :
Client Sample ID :
Date of Draw :
Date Received :
Reporting Date :
Resample :

Patient Information

Name :
Date of Birth :
Gest. Age at Draw :
Pregnancy type :
Indication :
Medical Record/Patient ID :

Provider Information

Hospital :
Physician :
Phone :
Fax :

Quality Test

Sample suitability	Pass	NGS data quality	Pass
DNA quality	Pass	Reference material test	Pass
Library quality	Pass	Fetal fraction	9.8%

Results

Chromosome	Result	PPV or NPV	Risk score (before)	Risk score (after)
Trisomy 21	Low Risk	>99(NPV)	1/307	<2/10,000
Trisomy 18	Low Risk	>99(NPV)	1/1,047	<2/10,000
Trisomy 13	Low Risk	>99(NPV)	1/3,222	<2/10,000
Trisomy 9	Low Risk	NA	NA	<2/10,000
Trisomy 16	Low Risk	NA	NA	<2/10,000
Trisomy 22	Low Risk	NA	NA	<2/10,000
XO	Low Risk	>99(NPV)	NA	<3/10,000
XXX	Low Risk	>99(NPV)	NA	<3/10,000
XXY	Low Risk	>99(NPV)	NA	<3/10,000
XXY	Low Risk	>99(NPV)	NA	<1/10,000
All other autosomal trisomies	Low Risk	NA	NA	<2/10,000

1. All probabilities and PPVs and NPVs are calculated on the site (<https://www.perinatalquality.org/Vendors/NSGC/NIPT/>) based on maternal age and NICE sensitivity/specificity data.
 • NPV : Negative Predictive Values (In case of Low Risk)
 • PPV : Positive Predictive Values (In case of High Risk)
2. The risk score (before) for aneuploidy was reported in a published study of 17,885 women [Dar et al. Am J Obstet Gynecol. 2014 Nov;211(5):527.e1-527.e17] based on maternal age, gestational age and/or general population.
3. The risk score (after) value is calculated as a probability value obtained through the Gaussian distribution of the normal group and the abnormal group using the fetal DNA percentage (fetal fraction), chromosome read% (Unique Read%) value, and Z-score value. The risk score (after) may not reflect the actual PPV of this patient, and other test results, ultrasound findings, and personal/family history are not included in the risk assessment.

Interpretation

Sample Report

NICE® Microdeletion/duplication Report



Sample Information

Sample Type :
Client Sample ID :
Date of Draw :
Date Received :
Reporting Date :
Resample :

Patient Information

Name :
Date of Birth :
Gest. Age at Draw :
Pregnancy type :
Indication :
Medical Record/Patient ID :

Provider Information

Hospital :
Physician :
Phone :
Fax :

Quality Test

Sample suitability	Pass	NGS data quality	Pass
DNA quality	Pass	Reference material test	Pass
Library quality	Pass		

Results

Location	Disease	Result	Location	Disease	Result
1p36	1p36 deletion syndrome	Low Risk	11q23	Jacobson syndrome	Low Risk
2q33.1	2q33.1 deletion syndrome	Low Risk	15q11.2-q13	Prader-willi / Angelman syndrome	Low Risk
4p16.3	Wolf-Hirschhorn syndrome	Low Risk	22q11.2	DiGeorge syndrome	Low Risk
5p-	Cri Du Chat syndrome	Low Risk	Etc	108 syndrome sites	High Risk
7q11.23	Williams-Beuren syndrome	Low Risk			

Interpretation

Sample Report

Limitations of Test

- This test is designed to screen for subchromosomal deletions in chromosomal regions: 1p36, 2q33.1, 4p16.3, 5p-, 7q11.23, 11q23, 15q11.2-q13, 22q11.2, 108 syndromes and is available for singleton pregnancies with gestational age of at least 10 weeks 0 days, as estimated by last menstrual period, crown rump length, or other appropriate method.
- These results do not eliminate the possibility that this pregnancy may be associated with other chromosomal or subchromosomal abnormalities, birth defects, and other conditions. This test is not intended to identify pregnancies at risk for open neural tube.
- In addition, there is a small possibility that the test results might not reflect the chromosome status of the fetus, but may reflect subchromosomal changes of the placenta (confined placental mosaicism), or of the mother.
- This is a screening test and this can result in false positive or false negative. Therefore negative results do not eliminate the possibility of 1p36 deletion, 2q33.1 deletion, 4p16.3 deletion, 5p- deletion, 7q11.23 deletion, 11q23 deletion, 15q11.2-q13 deletion, 22q11.2 deletion and 108 microdeletions/duplication syndromes. If definitive diagnosis is desired, chorionic villus sampling or amniocentesis would be necessary, with consideration of prenatal microarray or region specific DNA probes.
- In addition to the above-mentioned abnormalities, chromosomal and sub-chromosomal findings greater than 3 Mb may be reported. The advanced chromosomal and sub-chromosomal findings are rare and complex, insufficient validation may result in lower specificity. For microdeletion detection, Low fetal fraction or short CNV size is a technical limitation. Consultation on the advice of a physician(s) will be necessary for such findings.

List of 108 microdeletion/duplication syndromes

No.	Disease	No.	Disease	No.	Disease
1	1p32-q31 deletion syndrome	41	9p24.3 deletion syndrome	81	Miller-Dieker lissencephaly syndrome (MDLS) (loss)
2	1q41-q42 deletion syndrome	42	9q33.3q34.11 deletion syndrome	82	Miller-Dieker lissencephaly syndrome (MDLS) (gain)
3	1q43-q44 deletion syndrome	43	Early infantile epileptic encephalopathy 4 (EIEE4)	83	17p13.3 telomeric duplication syndrome
4	2p12-p11.2 deletion syndrome	44	Kleefstra syndrome 1 (KLEFS1)	84	17q21.31 deletion syndrome
5	2p15-p16.1 deletion syndrome	45	10p11.21-p12.31 microdeletion syndrome	85	17q21.31-q23.2 deletionsyndrome
6	2q13 deletion syndrome	46	DiGeorge syndrome/Velocardiofacial syndrome complex 2 (DSG2)	86	17q21.31-q23.2 deletionsyndrome
7	2q13 duplication syndrome	47	10q22.3-q23.2 deletion syndrome	87	Tetrasomy 18p syndrome
8	2q31.1 microdeletion syndrome	48	Split hand/foot malformation 3 (SHFM3)	88	18p deletion syndrome
9	2q31.1 duplication syndrome	49	10q26 deletion syndrome	89	18q deletion syndrome
10	2q35 duplication syndrome	50	Potocki-Shaffer syndrome	90	19p13 duplication syndrome
11	3p25.3 deletion syndrome	51	WAGR syndrome	91	19q13.11 microdeletion syndrome
12	3pter-p25 deletion syndrome	52	WAGRO syndrome	92	20p13 microdeletion syndrome
13	3q13.31 deletion syndrome	53	11q13.2-q13.4 deletion syndrome	93	21q22.11-q22.12 microdeletion syndrome
14	Dandy-Walker syndrome (DWS)	54	11q22.2-q22.3 microdeletion syndrome	94	22q11.2 deletion syndrome (distal, D-E/F)
15	3q26 microduplication syndrome	55	11q23 deletion syndrome	95	22q11.2 deletion syndrome (LCR22 B/C/D)
16	3q29 deletion syndrome	56	12p21.1 microdeletion syndrome	96	22q13 deletion syndrome
17	4q21 deletion syndrome	57	12q14 microdeletion syndrome	97	22q13 duplication syndrome
18	Aventfield-Rieger syndrome, type 1 (RIECS)	58	12q15q11.1 microdeletion syndrome	98	Xp11.22 duplication syndrome
19	5p13 duplication syndrome	59	13q14 deletion syndrome	99	Xp11.22-p11.23 duplication syndrome
20	5q12 deletion syndrome	60	14q11-q22 deletion syndrome	100	Xp11.23 microdeletion syndrome
21	5q14.3 deletion (proximal) syndrome	61	Frias syndrome	101	Xp11.3 deletion syndrome
22	Sotos syndrome	62	14q24.1-q24.3 microdeletion syndrome	102	Xp21 microdeletion syndrome
23	6p22 microdeletion syndrome	63	15q13.3 deletion syndrome (BP4 to BPS) (loss)	103	Xp21 microduplication syndrome
24	6q11-q14 deletion syndrome	64	15q13.3 deletion syndrome (BP4 to BPS) (gain)	104	Xp22.31 microdeletion syndrome
25	6q24-q25 deletion syndrome	65	15q14 microdeletion syndrome	105	Xq21 microdeletion syndrome
26	Coffin-Siris syndrome 1 (CSS1)	66	15q25.2 deletion (proximal) syndrome	106	Xq22.3 telomeric deletion syndrome
27	Chordoma	67	15q26-qter deletion syndrome	107	Xq27.3-q28 duplication syndrome
28	Greig cephalopolysyndactyly syndrome (GCPSS)	68	16p11.2-p12.2 microduplication syndrome	108	Xq28 deletion syndrome
29	7p22.1 microduplication syndrome	69	16p12.2 deletion (proximal) syndrome		
30	7q11.23 deletion (distal) syndrome	70	16p13.11 duplication syndrome		
31	Williams-Beuren syndrome (WBS)	71	16p13.11 deletion syndrome		
32	Curranio syndrome	72	Polycystic kidney disease, infantile severe, with tubular sclerosis (PKDTS)		
33	7q36.3 duplication syndrome	73	Rubinstein-Taybi syndrome		
34	8p11.2 deletion syndrome	74	Alpha-thalassemia/intellectual disability syndrome, chromosome 16 related (ATR 16 syndrome)		
35	8p23.1 deletion syndrome	75	16q22 duplication syndrome		
36	8q12 microduplication syndrome	76	Smith-Magenis syndrome		
37	Nablas mask-like facial syndrome (NMLFS)	77	Yuan-Hareé-Lupski syndrome (YUHAL)		
38	Trichorhinopharyngeal syndrome type 2 (TRPS2)	78	17p12 deletion syndrome		
39	9p deletion syndrome	79	17p12 duplication syndrome		
40	9p13 microdeletion syndrome	80	17p13.1 deletion syndrome		



**Available after
10 weeks pregnant**



**More than 99%
test success rate**






**Safe and easy test
with maternal blood**






**Bioinformatics
pipeline accuracy**

Chromosomal heatmap






10mb bin heatmap






Original

Fetus




T9 17W 6D
FF : 5.5
47,X(),+9[14]/46,X()[16]



T16 17W 3D
FF : 14.98
47,X(),+16[3]/46,X()[37]




T22 17W 0D
FF : 10.92
47,X(),+22



* Total : **28,800** samples (Artificial samples due to lack of clinical samples)
 * Create a known region for DECIPHER and OMIM

Disease	deletion location and length (DECIPHER)	Detected Range (Mb)	LOD length (Mb)	LOD Fetal Fraction	sensitivity	specificity
1p36 deletion syndrome	1:10001-12840259 (12.83Mb)	0.5~12.83	≥3	≥3%	42~100%	100%
2q33.1 deletion syndrome	2:196925121-205206939 (8.23Mb)	0.5~10	≥3	≥3%	20~97%	100%
Wolf-Hirschhorn syndrome (4p16.3 deletion)	4:1569197-2110236 (0.54Mb)	1~10	≥3	≥5%	19~87%	100%
Cri du chat syndrome (5p15.3 deletion)	5:10001-12533304 (12.52Mb)	0.5~12.52	≥4	≥5%	44~95%	100%
Williams beuren syndrome (7q11.23 deletion)	7:72744455-74142672 (1.39Mb)	0.5~10	≥4	≥3%	46~92%	100%
Jacobsen syndrome (11q23 deletion)	11:110470724-121170709 (10.69Mb)	0.5~10.69	≥3	≥3%	27~96%	100%
Prader-Willi / Angelman syndrome (15q11.2 deletion)	15:22749354-28438266 (5.68Mb)	0.5~10	≥4	≥3%	20~87%	100%
DiGeorge syndrome (22q11.2 deletion)	22:19009792-21452445 (2.44Mb)	2~10	≥5	≥5%	34~97%	100%



- Technical limitation : Low Fetal Fraction, Short CNV size
- In the shallow-depth NGS sequencing method, there are chromosomal regions that are difficult to read mapping to the reference sequence depending on the chromosomal characteristics.